

the Fungitell® Bulletin

volume 10, issue 3

Topic:

CANDIDA AURIS AND (1→3)-β-D-GLUCAN


Fungitell® Bulletins are intended as technical advisory communications and as such are disseminated to the general public in order to highlight the significance of (1→3)-β-D-Glucan on human health. These communications do not promote a specific drug, therapy nor make any representation or suggestion concerning the suitability or effectiveness of a particular drug or therapy in patients harboring (1→3)-β-D-Glucan. Fungitell® is an adjunct diagnostic assay to be utilized in conjunction with clinical signs and symptoms for the diagnosis of invasive fungal infection. Fungitell® is currently 510(k) cleared for the detection and quantification of (1→3)-β-D-Glucan in human serum and should be used and interpreted only in a manner consistent with the current Instructions for Use.

Discussion:

The (1→3)-β-D-glucan (BDG) assay is widely recognized for its high negative predictive value (NPV)^{1,2}. When used in conjunction with additional evidence, the BDG assay allows for the exclusion of infections caused by various fungal organisms, including *Candida*, *Aspergillus*, *Pneumocystis*, and most endemic fungi. A negative BDG test in high-risk patients with suspected invasive fungal disease provides clinicians with confidence to modify empiric antifungal therapy, leading to improved patient outcomes³. However, when positive BDG results are reported, it becomes important to understand the varying sensitivity of the assay, especially in cases where certain species of microorganisms are highly prevalent. This may be relevant in the context of the significant increase in global cases of *Candida auris*.

Candida auris is a multidrug-resistant yeast that can cause severe invasive infections, particularly in hospitalized patients with underlying health conditions. Classified as a 'critical priority'⁴ and an 'urgent threat'⁵, *Candida auris* has garnered increased attention due to a 61% rise in clinical cases in the US from 1,474 (2021) to 2,377 (2022)⁶. First isolated almost 15 years ago⁷, cases of *Candida auris* have now been reported in thirty-nine countries⁸. The mortality rate for *Candida auris* infections is high, ranging from 30-72%, with the most severe invasive infections occurring in critically ill patients⁹. Although several meta-analyses have investigated the use of Fungitell for the detection of invasive fungal disease in a critical care setting¹⁰, there are limited publications specifically examining the utility of Fungitell in relation to *Candida auris* infection.

Timely identification of patients is crucial for effective *Candida auris* infection prevention and control¹¹. Parak *et al.*¹² demonstrated that all cases of *Candida auris* infection in their study were healthcare-acquired (positive culture obtained >48 hours after admission). Chibabhai *et al.*¹³ showed that a positive BDG (≥80pg/mL) was observed approximately 2.5 days **prior** to blood culture collection, showing utility of timely serial BDG measurements in monitoring high-risk patients for candidemia. This finding is supported by Tisot



Bulletin Volume 10, issue 3
Publish Date: August 2023
CORP_0516

Corporate Headquarters
Associates of Cape Cod, Inc.
124 Bernard E. Saint Jean Drive
East Falmouth, MA 02536 USA
Tel: (508) 540-3444
www.acciusa.com

United Kingdom
Associates of Cape Cod Int'l., Inc.
Deacon Park, Moorgate Road
Knowsley, Liverpool L33 7RX
United Kingdom
Tel: (44) 151-547-7444
www.acciu.co.uk

Europe
Associates of Cape Cod Europe GmbH
Opelstrasse 14
D-64546 Mörfelden-Walldorf, Germany
Tel: (49) 61 05-96 10 0
www.acciusa.de

et al.¹⁴ which showed early BDG positivity in relation to other *Candida* species.

Levels of serum BDG in confirmed cases of *Candida auris* were found to be lower (129pg/mL) compared to *Candida albicans* (410pg/mL)^{15,16}. The reduced BDG burden in *Candida auris* candidemia leads to decreased BDG sensitivity (from ~80% to 71-43%) compared to candidemia caused by other *Candida* species^{12,13,16}. Affinity for biofilm formation, cell wall composition among *Candida* species, and prior antifungal therapy, are believed to influence BDG sensitivity. Furthermore, detectable BDG levels may vary among different clades of *Candida auris*¹⁶, illustrated by the variability in reported serum BDG levels in studies from South Africa¹³, Pakistan¹⁷, and Italy¹⁶.

In summary, lower BDG values may be observed when screening patients in intensive care with a high prevalence of *Candida auris* infection. In cases of suspected candidemia, measuring BDG levels at enrolment and subsequently every 48-72 hours for 14 days is recommended to increase sensitivity. The high negative predictive value of Fungitell can guide antimicrobial therapy, leading to a reduction in the duration of empirical antifungal treatment for critically ill septic/candidemic patients¹⁸. Despite the decreased sensitivity observed in some patients, elevated BDG levels may still serve as a useful marker for *Candida auris* candidemia.

Discussion References

- Donato, L., González, T., Canales, M., Legarraga, P., García, P., & Rabagliati, R. (2017). Evaluación del rendimiento de 1,3-β-D-glucano como apoyo diagnóstico de infecciones invasoras por *Candida* spp. en pacientes críticos adultos [The 1,3-β-D-glucan in critical adult patients as diagnostic tool for invasive *Candida* spp. infection, performance evaluation]. *Revista chilena de infectología: organo oficial de la Sociedad Chilena de Infectología*, 34(4), 340–346. <https://doi.org/10.4067/s0716-10182017000400340>
- Giacobbe, D. R., Mikulska, M., Tumbarello, M., Furfaro, E., Spadaro, M., Losito, A. R., Mesini, A., De Pascale, G., Marchese, A., Bruzzone, M., Pelosi, P., Mussap, M., Molin, A., Antonelli, M., Posteraro, B., Sanguinetti, M., Viscoli, C., Del Bono, V., & ISGRI-SITA (Italian Study Group on Resistant Infections of the Società Italiana Terapia Antinfettiva) (2017). Combined use of serum (1,3)-β-D-glucan and procalcitonin for the early differential diagnosis between candidaemia and bacteraemia in intensive care units. *Critical care (London, England)*, 21(1), 176. <https://doi.org/10.1186/s13054-017-1763-5>
- Rautemaa-Richardson, R., Rautemaa, V., Al-Wathiqi, F., Moore, C. B., Craig, L., Felton, T. W., & Muldoon, E. G. (2018). Impact of a diagnostics-driven antifungal stewardship programme in a UK tertiary referral teaching hospital. *The Journal of antimicrobial chemotherapy*, 73(12), 3488–3495. <https://doi.org/10.1093/jac/dky360>
- WHO. WHO fungal priority pathogens list to guide research, development and public health action. Geneva: World Health Organization; 2022. 2022.
- CDC. CDC. 2019. Antibiotic Resistance Threats in the United States, 2019. Atlanta, GA: U.S. Department of Health and Human Services, CDC; 2019.
- CDC. <https://www.cdc.gov/fungal/candida-auris/tracking-c-auris.html>. [Online] 14 February 2023.
- Satoh, K., Makimura, K., Hasumi, Y., Nishiyama, Y., Uchida, K., & Yamaguchi, H. (2009). *Candida auris* sp. nov., a novel ascomycetous yeast isolated from the external ear canal of an inpatient in a Japanese hospital. *Microbiology and immunology*, 53(1), 41–44. <https://doi.org/10.1111/j.1348-0421.2008.00083.x>
- Sanyaolu, A., Okorie, C., Marinkovic, A., Abbasi, A. F., Prakash, S., Mangat, J., Hosen, Z., Haider, N., & Chan, J. (2022). *Candida auris*: An Overview of the Emerging Drug-Resistant Fungal Infection. *Infection & chemotherapy*, 54(2), 236–246. <https://doi.org/10.3947/ic.2022.0008>
- Cortegiani, A., Misseri, G., Fasciana, T., Giammanco, A., Giarratano, A., & Chowdhary, A. (2018). Epidemiology, clinical characteristics, resistance, and treatment of infections by *Candida auris*. *Journal of intensive care*, 6, 69. <https://doi.org/10.1186/s40560-018-0342-4>
- Xiaoling, L., Tingyu, T., Caibao, H., Tian, Z., & Changqin, C. (2018). Diagnostic Efficacy of Serum 1,3-β-D-glucan for Invasive Fungal Infection: An Update Meta-Analysis Based on 37 Case Or Cohort Studies. *Open medicine (Warsaw, Poland)*, 13, 329–337. <https://doi.org/10.1515/med-2018-0050>
- Lyman, M., Forsberg, K., Sexton, D. J., Chow, N. A., Lockhart, S. R., Jackson, B. R., & Chiller, T. (2023). Worsening Spread of *Candida auris* in the United States, 2019 to 2021. *Annals of internal medicine*, 176(4), 489–495. <https://doi.org/10.7326/M22-3469>
- Parak, A., Stacey, S. L., & Chibabhai, V. (2022). Clinical and laboratory features of patients with *Candida auris* cultures, compared to other *Candida*, at a South African Hospital. *Journal of infection in developing countries*, 16(1), 213–221. <https://doi.org/10.3855/jidc.14917>
- Chibabhai, V., Fadana, V., Bosman, N., & Nana, T. (2019). Comparative sensitivity of 1,3 beta-D-glucan for common causes of candidaemia in South Africa. *Mycoses*, 62(11), 1023–1028. <https://doi.org/10.1111/myc.12982>
- Tissot, F., Lamoth, F., Hauser, P. M., Orasch, C., Flückiger, U., Siegemund, M., Zimmerli, S., Calandra, T., Bille, J., Eggmann, P., Marchetti, O., & Fungal Infection Network of Switzerland (FUNGINOS) (2013). β-glucan antigenemia anticipates diagnosis of blood culture-negative intraabdominal candidiasis. *American journal of respiratory and critical care medicine*, 188(9), 1100–1109. <https://doi.org/10.1164/rccm.201211-2069OC>
- Mikulska, M., Giacobbe, D. R., Furfaro, E., Mesini, A., Marchese, A., Del Bono, V., & Viscoli, C. (2016). Lower sensitivity of serum (1,3)-β-D-glucan for the diagnosis of candidaemia due to *Candida parapsilosis*. *Clinical microbiology and infection: the official publication of the European Society of Clinical Microbiology and Infectious Diseases*, 22(7), 646.e5–646.e6.46E8. <https://doi.org/10.1016/j.cmi.2016.05.020>
- Mikulska, M., Furfaro, E., Magnasco, L., Coddà, G., Giacobbe, D. R., Dentone, C., Vena, A., Marchese, A., & Bassetti, M. (2022). Levels of beta-D-glucan in *Candida auris* supernatants, an in vitro and in vivo preliminary study. *Clinical microbiology and infection: the official publication of the European Society of Clinical Microbiology and Infectious Diseases*, 28(8), 1154.e1–1154.e3. <https://doi.org/10.1016/j.cmi.2022.02.045>
- Farooqi, J., Niamatullah, H., Irfan, S., Zafar, A., Malik, F., & Jabeen, K. (2021). Comparison of β-D-Glucan levels between *Candida auris* and other *Candida* species at the time of candidaemia: a retrospective study. *Clinical microbiology and infection: the official publication of the European Society of Clinical Microbiology and Infectious Diseases*, 27(10), 1519.e1–1519.e5. <https://doi.org/10.1016/j.cmi.2021.05.031>
- De Pascale, G., Posteraro, B., D'Arrigo, S., Spinazzola, G., Gaspari, R., Bello, G., Montini, L. M., Cutuli, S. L., Grieco, D. L., Di Gravio, V., De Angelis, G., Torelli, R., De Carolis, E., Tumbarello, M., Sanguinetti, M., & Antonelli, M. (2020). (1,3)-β-D-Glucan-based empirical antifungal interruption in suspected invasive candidiasis: a randomized trial. *Critical care (London, England)*, 24(1), 550. <https://doi.org/10.1186/s13054-020-03265-y>

